

# URICHECK

## 11 parameters

REF U003C

150


Reagent strips for urinalysis  
Bandelettes reactives pour analyse d'urine  
Tiras reactivas para analisis de orina  
Reagenz streifen für Urinanalyse

CE

Replace cap promptly and tightly.  
Refermer rapidement et soigneusement le flacon après usage.  
Cerrar el bote rápidamente.

IVD

Dose nach jeder Entnahme sofort wieder fest verschließen.



VEDA LAB

ZAT du Londeau  
Rue de l'Expansion - Cerisé  
BP 181 - 61006 ALENCON CEDEX - FRANCE

Tests	Results / Resultats / Resultados / Ergebnisse
Blood/Sang Sangre/Blut	
Bilirubin/Bilirubine Bilirubina	
Urobilinogen/Urobilinogène Urobilinógeno	
Ketones/Cétones Cetonas/Ketonkörper	
Glucose/Glucosa Glukose	
Protein/Protéines Proteínas	
Nitrite/Nitritos Nitrit	
Leukocytes/Leucocytes Leucocitos/Leukozyten	
pH	
S.G./Densité Densidad/Spesz.Gew.	
Ascorbic acid/Acide ascorbique Acido Ascórbico/Ascorbinsäure	

LOT



Read at exactly 60 sec. (Leukocytes 90-120 sec.)  
Lire après 60 sec. (Leucocytes 90-120 sec.)  
Leer exactamente a los 60 seg. (Leucocitos 90-120 seg.)  
Nach exakt 60 Sek. ablesen (Leukozyten 90-120 sek.)